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1. Immunoprecipitation

The tissue was washed with PBS, followed by the addition of lysis buffer in three times the volume of the tissue. The tissue was then homogenized using bead milling, and subsequently subjected to sonication to disrupt the cells. The mixture was incubated on ice for 30 min, followed by centrifugation to collect the supernatant. The supernatant was transferred to a new EP tube, and the protein concentration was determined.. Vortex magnetic bead stock solution, transfer bead suspension to EP tube, remove supernatant with MPC, wash beads with 1 mL PBST, and resuspend beads in 500 μ L PBS. Add sufficient antibody for immunoprecipitation and rotate at room temperature (18-22° C) for 1 hour. Next, remove supernatant using MPC, wash four times with 1 mL PBS. Add antibody-coated beads to the cell lysate sample, incubate overnight at 4° C with rotation, then wash thoroughly with PBST. For the final wash, transfer beads to a fresh tube, remove supernatant, and elute proteins with protein loading buffer. Run eluted samples on SDS-PAGE. Transfer one gel for subsequent Western blot analysis. Stain another gel with R250 Coomassie Brilliant Blue for 0.5 hours, then perform gel strip detection.

2. In-gel Digestion

For in-gel tryptic digestion, gel pieces were destained in 50 mM NH_4HCO_3 in 50% acetonitrile (v/v) until clear. Gel pieces were dehydrated with 100 μl of 100% acetonitrile for 5 min, the liquid removed, and the gel pieces rehydrated in 10 mM dithiothreitol and incubated at 56 °C for 60 min. Gel pieces were again dehydrated in 100% acetonitrile, liquid was removed and gel pieces were rehydrated with 50 mM iodoacetamide. Samples were incubated at room temperature, in the dark for 45 min. Gel pieces were washed with 50 mM NH_4HCO_3 and dehydrated with 100% acetonitrile. Gel pieces were rehydrated with 10 ng/ μl trypsin resuspended in 50 mM NH_4HCO_3 and 5 mM CaCl_2 on ice for 10 min. Excess liquid was removed and gel pieces were digested with trypsin at 37 °C overnight. Peptides were extracted with 50% acetonitrile /5% formic acid, followed by 100% acetonitrile. Peptides were dried to completion and resuspended in 2% acetonitrile/0.1% formic acid.

3. LC-MS/MS Analysis

The tryptic peptides were dissolved in 0.1% formic acid (solvent A), directly loaded onto a home-made reversed-phase analytical column (15-cm length, 75 μm i.d.). The gradient was comprised of an increase from 5% to 34% solvent B (0.1% formic acid in 80% acetonitrile) over 40 min, 34% to 38% in 5 min , 38% to 90% in 10 min and then holding at 90% for the last 10 min, all at a constant flow rate of 300 nl/min on an EASY-nLC 1200 UPLC system.

The peptides were subjected to NSI source followed by tandem mass spectrometry (MS/MS) in Thermo Scientific Orbitrap Exploris 480 coupled online to the UPLC. The electrospray voltage applied was 2.0 kV. The m/z scan range was 350 to 1800 for full scan, and intact peptides were detected in the Orbitrap at a resolution of 70,000. Peptides were then selected for MS/MS using NCE setting as 28 and the fragments

were detected in the Orbitrap at a resolution of 17,500. A data-dependent procedure that alternated between one MS scan followed by 20 MS/MS scans with 15.0s dynamic exclusion. Automatic gain control (AGC) was set at 5E4.

3. Data Processing

The resulting MS/MS data were processed using Proteome Discoverer 2.4. Tandem mass spectra were searched against *Uniprot _Mus musculus*. Trypsin/P was specified as cleavage enzyme. Mass error was set to 10 ppm for precursor ions and 0.02 Da for fragment ions.